

VLIFT® Implant Sizing

DESIGNATION	HEIGHT RANGE					Distraction Range
	Undistracted No End Caps	Fully Distracted No End Caps	Undistracted Two, 0° End Caps	Fully Distracted Two, 0° End Caps		
Ø18mm x 20mm	18.5	25.5	20.5	27.5	7	
Ø18mm x 25mm	23	34.5	25	36.5	11.5	
Ø18mm x 32mm	30	48.5	32	50.5	18.5	
Ø18mm x 32mm + 1 extension	45	63.5	47	65.5	18.5	
Ø18mm x 32mm + 2 extensions	60	78.5	62	80.5	18.5	
Ø22mm x 25mm	23	34.5	25	36.5	11.5	
Ø22mm x 32mm	30	48.5	32	50.5	18.5	
Ø22mm x 37mm	35	58.5	37	60.5	23.5	
Ø22mm x 37mm + 1 extension	50	73.5	52	75.5	23.5	
Ø22mm x 37mm + 2 extensions	65	88.5	67	90.5	23.5	

Note: Dimensions are represented at the undistracted and fully distracted mean values.

All dimensions are in (mm).

